

Visually Lossless Compression of Chest X-ray Images

Mohammad Bagher Sharifi
Ali Naseri Esfahani
Erfan Moeini



Overview

- **Introduction:** prior work & project overview
- **Method:** noise analysis & compression approach
- **Results:** dataset, experiments, outcomes

Introduction

Related Work

A Fast Method of Visually Lossless Compression of Dental Images


- Dental images
- Visually-lossless
- Compression efficiency



Goal


Determine optimal Quantization Step (QS)

- Ensure visually-lossless output
- Maximize compression ratio without harm




Framework

- Noise analysis and modeling
- Coder comparison using perceptual metrics
- QS selection strategy
- Clinical verification



Our Work

- Implement article's compression
- Reproduce experiments
- Test on chest X-rays



Method




Methodology

- This work applies the ADCT algorithm to compress images
- Because ADCT adapts compression to image content, it provides better fidelity than non-adaptive methods such as JPEG2000
- This work aims to select the optimal input hyperparameter for the ADCT algorithm

ADCT Algorithm


- based on DCT with adaptive partitioning
- driven by Quantization Step (QS) as a key input parameter
- context-based statistical modeling and arithmetic coding of DCT data
- choosing QS balances compression ratio and diagnostic quality



Analysis of Noise Properties in Dental Images

- Noise variance is estimated as a combination of additive, multiplicative, and Poisson-like components

$$\hat{\sigma}_{eq}^2 = \hat{\sigma}_a^2 + \hat{\sigma}_\mu^2 I^2 + kI$$




Invisibility Condition of Compression Noise in Images

- For compression distortion to be invisible, it must be much smaller than the image noise

$$\sigma_{distortion}^2 \leq \frac{\sigma_{noise}^2}{10} \quad (1)$$

$$PSNR = 10 \log_{10} \left(\frac{MAX^2}{MSE^2} \right) \quad (2)$$





Methodology Overview

Goal: Efficient compression of dental images with high diagnostic quality

Idea: Select the optimal Quantization Step (QS) based on image noise

Workflow:

1. Noise analysis of input image
2. Define target PSNR from noise level
3. Apply ADCT with optimal QS
4. Evaluate quality with metrics

Result



Experiments

- **Datasets:** panoramic dental images and chest images in two resolutions
- **Metrics:** FSIM, PSNR HVS M
- **Experiments:**
 - Validate paper results on a new dental dataset
 - Find QS for chest X-ray Images
 - Analyze QS in higher resolution

Table1. Datasets Details

Type	Imaging	Number Of Images
Panoramic Dental [1]	X-Ray 512	48
Chest [2]	X-Ray 512	145
Chest [2]	X-Ray 1024	107

ADCT vs JPEG2000

- Encode and decode Images with different QS values
- Encoding using jpeg2000
- Significant improvement compared to jpeg2000



Fig1. Sample of panoramic dental image




Fig2. adct and jpeg2000 quality on dental image₁₆

Find QS

- Encode and decode images for QS in the range (6, 42)
- Find QS values where
 - FSIM > 0.99
 - Psnr hvsm > 40db
- Charts
 - Red line shows the threshold
 - Blue line indicates the best QS for each metric
- QS=16 works well for all sample images




Fig3. Adct on panoramic dental dataset

Best QS for Chest X-ray 512

- QS = 14 performs well in psnr-hvs-m
- QS = 6 performs well in fsim
- The average CR is around 10
- Results are not always normal




Fig4. image with low fsim




Fig5. Adct on chest X-ray 512 dataset 18

Best QS for Chest X-ray 1024

- QS = 12 also performs well in psnr-hvs-m
- QS = 6 does not perform well in fsim
- A single good QS value could not be found for all scenarios




Fig6. Adct on chest X-ray 1024 dataset 19



Clinical Verification

- Chest X-ray 512 compress with QS=6
- Image shuffling and doctor evaluation
- Better in Ideal images

Table2. results of doctor evaluation

Compression Status	Total Images	Unusable %	Weak %	Sufficient %	Ideal %
Compressed	27	0	22	33	45
Uncompressed	23	4	22	35	39



Conclusion

- **Introduction:** Main paper overview
- **Method:** Relation between noise and qs , Explain Adct
- **Results:** Check paper results and Find QS for chest images

Any Question?



References

1. Mooney, P. (2018). Chest X-Ray Images (Pneumonia) [Data set]. Kaggle. Retrieved August 19, 2025, from <https://www.kaggle.com/datasets/paultimothymooney/chest-xray-pneumonia>
2. Vipul Arya, K. (n.d.). A Collection of Dental X-Ray Images for Analysis [Data set]. Kaggle. Retrieved August 19, 2025, from <https://www.kaggle.com/datasets/kvipularya/a-collection-of-dental-x-ray-images-for-analysis>
3. Krivenko, S., Lukin, V., Krylova, O., Kryvenko, L., & Egiazarian, K. (2021). A Fast Method of Visually Lossless Compression of Dental Images. **Applied Sciences**, 11(1), Article 135. <https://doi.org/10.3390/app11010135>
4. Abramova, V., Krivenko, S., Lukin, V., & Krylova, O. (2020). Analysis of noise properties in dental images. In **2020 IEEE 40th International Conference on Electronics and Nanotechnology (ELNANO)** (pp. 511–515). IEEE. <https://doi.org/10.1109/ELNANO50318.2020.9088819>
5. Ponomarenko, N., Lukin, V., Egiazarian, K., & Astola, J. (2008). ADCTC: Advanced DCT-Based Image Coder. In **Proceedings of the 2008 International Workshop on Local and Non-Local Approximation in Image Processing (LNLA 2008)** (pp. 511–515). IEEE. <https://doi.org/10.1109/ELNANO50318.2008.4522211>

Appendix



Datasets

Table1. Datasets Details

Type	Imaging	Resolution	Number Of Images	Sample
Panoramic Dental	X-Ray	512 * 512	48	35
Chest	X-Ray	512 * 512	145	26
Chest	X-Ray	1024 * 1024	107	70